# Zhlus®

# ATLANTIC WORLD "Bringing Security to life"

### MEDICA GAUZE SWAB B.P. (Sterile)

### SIZE:

5 CM X 5 CM X 16 PLY, 7.5 CM X 7.5 CM X 16 PLY, 10 CM X 10 CM X 16 PLY

### CHARACTERISTICS

- · Made From Accomplished Cotton Yarn
- · Absorption rate faster (less than 10 seconds)
- · Better tending capacity
- · Less allegiance to healing tissue
- · These Gauze Swab available in American fold and British fold
- · Sterilized by gamma radiation

### **APPLICATION**

- It is an useful utility product in kit, can be used for padding & protection
- · For swabbing and cleaning wounds
- · For dressing low exudates wounds



### Medica X-Ray Detectable Abdominal Pad B.P.

### SIZE:

25 CM X 25 CM X 8 PLY, 25 CM X 40 CM X 8 PLY, 40 CM X 40 CM X 6 PLY

### **CHARACTERISTICS**

- Made From Accomplished 100 % Cotton Yarn
- Faster rate of absorption (Less than 10 Sec)
- · Higher tending capacity · Sterilized by Gamma
- Sponge contains: Layers of highly absorbent cotton gauze stitched together with X-ray detectable thread and sturdy loop attached on the one side of corner
- · X-Ray Thread Help in Detect the Presence of Unwanted Things
- · Also available in Sterile EO

- · During Surgery used to absorb blood & exudates
- · Clean & prevent wounds from being infected





# Zplus®



### EYE PAD B.P.

### Medica Eye Pad

SIZE:

STD (6 CM X 8 CM)

### **CHARACTERISTICS**

- Highly absorbent
- Extra soft and aseptic
- · Protects from air, light & dust
- Sterilized by gamma radiation

### **APPLICATION**

For eye dressing



### **MEDIPAD**

### Medica Combine (operation) Surgical Pad

SIZE:

10 CM x 10 CM, 10 CM x 20 CM, 10 CM x 15 CM

### **CHARACTERISTICS**

- Softer
- Dressing pad is made from absorbent cotton + non-woven cloth
- Faster rate of absorption and more tending capacity
- Sterilized by gamma radiation
- Also available in Sterile EO

### **APPLICATION**

- Better Care & Supportive Role in Wound Dressing & Surgery
- For aseptic post-operative dressings
- · Keep the Plain side on the operated region / wound and affix the adhesive plaster



### **MEDIFIN**

### Medica Paraffin Gauze Dressing B.P.

### SIZE:

10 CM x 10 CM X 01 NOS. 10 CM x 10 CM X 02 NOS. 10 CM x 30 CM X 01 NOS. 10 CM x 10 CM X 05 NOS. 10 CM x 10 CM X 10 NOS. 10 CM x 40 CM X 01 NOS. 10 CM x 10 CM X 20 NOS. 10 CM x 10 CM X 24 NOS. 10 CM x 10 CM X 36 NOS. 10 CM x 07 MTR X 01 NOS.

### **CHARACTERISTICS**

- Paraffin acts as a soothing agent. The antibiotic / antiseptic in it helps in the healing process.
- Interlocked gauze yarn threads, minimize when dressing is cut to shape.
- Does not stick to the wound and does not disturb the dressing process.
- Due to special type of butter paper loss of medicine after removal of butter paper is minimum.
- Sterilized by gamma radiation

- On burnt skin and prior to any dressings
- For scalds, leg ulcers and on skin graft sites, plastic & cosmetic surgery
- Reduces risk of colonization







### **CHLODINE**

Chlorhexidine Gauze Dressing B.P.
Contain 0.5% W/W Chlorhexidine Acetate B.P

### SIZE:

10 CM x 10 CM, 10 CM x 30 CM, 10 CM x 40 CM

### **CHARACTERISTICS**

- It contains 0.5% W/W chlorhexidine Acetate B.P.
- Chlorhexidine Gauze Dressing acts as a soothing agent. The antibiotic / antiseptic in it helps in the healing process.
- Interlocked gauze yarn threads, minimize when dressing is cut to shape.
- Due to special type of butter paper, loss of medicine after removal of butter paper is minimum.
- Sterilized by gamma radiation

### APPLICATION

- · On burnt skin and prior to any dressings
- For scalds, leg ulcers and on skin graft sites, plastic & cosmetic surgery
- Reduces risk of colonization



Gauze Swab B.P. Sterilized by EO

### SIZE:

5 CM x 5 CM x 16 PLY X 1 NOS. (20 CM x 20 CM UNFOLDED)
7.5 CM x 7.5 CM x 16 PLY X 1 NOS. (30 CM x 30 CM UNFOLDED)
10 CM x 10 CM x 16 PLY X 1 NOS (40 CM x 40 CM UNFOLDED)
5 CM X 5 CM FOLDED IN 8 PLY/12 PLY, 4 Pieces in Pack,
7.5 CM X 7.5 CM, FOLDED IN 8 PLY/12 PLY, 4 Pieces in Pack,
10 CM X 10 CM, FOLDED IN 8 PLY/12 PLY, 4 Pieces in Pack

### **CHARACTERISTICS**

- Made from Accomplished Cotton Yarn
- Absorption rate faster (less than 10 seconds)
- · Better tending capacity
- · Less allegiance to healing tissue
- · These Gauze Swab available in American fold and British fold
- Also available with X-ray detectable thread

- It is an useful utility product in kit, can be used for padding & protection
- · For swabbing and cleaning wounds
- · For dressing low exudates wounds









### **MEDIRIBBON**

### Medica Absorbent Cotton Ribbon Gauze B.P.

### SI7F:

2.5 CM X 2 MTR, 2.5 CM X 5 MTR

### CHARACTERISTICS

- Made from Accomplished 100 % cotton Yarn
- · Higher Rate of Absorption
- · Folded in four fold
- · Sterile By gamma Radiation / EO

### **APPLICATIONS**

- · For Gynec Surgery
- E.N.T. Dressing
- · Intended for use in management of clean
- Used in contaminated exuding wound like cavity wound & fistulas



### **MEDIMOP**

### X -Ray Detectable Abdominal Pad B.P.

### SI7F:

10 CM x 10 CM x 8 PLY, 15 CM x 15 CM x 8 PLY, 25 CM x 25 CM x 8 PLY, 30 CM x 30 CM x 8 PLY, 45 CM x 45 CM x 8 PLY, 37 CM x 45 CM x 8 PLY, 10 CM x 10 CM x 12 PLY, 15 CM x 15 CM x 12 PLY, 25 CM x 25 CM x 12 PLY, 30 CM x 30 CM x 12 PLY, 45 CM x 45 CM x 12 PLY, 37 CM x 45 CM x 12 PLY, 37 CM x 45 CM x 6 PLY

### **CHARACTERISTICS**

- Made From Accomplished 100 % Cotton Yarn
- Faster rate of absorption (Less than 10 Sec)
- · Higher tending capacity
- Sponge contains: Layers of highly absorbent cotton gauze stitched together with X-ray detectable thread and sturdy loop attached on the one side of corner
- X-Ray Thread Help in Detect the Presence of Unwanted Things

### **APPLICATIONS**

- · During Surgery used to absorb blood & exudates
- Clean wounds & prevent wounds from being infected

# MEDIRIBON

### Absorbent Cotton Ribbon Gauze B.P.

### SIZE:

1.25 CM X 100 MTR, 2.50 CM X 100 MTR

### **CHARACTERISTICS**

- Made from Accomplished 100 % cotton Yarn
- · Higher Rate of Absorption
- · Folded in four fold

- For Gynec Surgery
- · E.N.T. Dressing
- · Intended for use in management of clean
- Used in contaminated exuding wound like cavity wound & fistulas









### **GAUZE SWAB B.P.**

### Medica Gauze Swab B.P.

### SIZE:

5 CM x 5 CM x 8 PLY, 7.5 CM x 7.5 CM x 8 PLY, 10 CM x 10 CM x 8 PLY, 5 CM x 5 CM x 12 PLY, 7.5 CM x 7.5 CM x 12 PLY, 10 CM x 10 CM x 12 PLY, 7.5 CM x 7.5 CM x 12 PLY, 10 CM x 10 CM x 12 PLY

### **CHARACTERISTICS**

- · Made From Accomplished Cotton Yarn
- · Absorption rate faster (less than 10 seconds)
- · Better tending capacity
- · Less allegiance to healing tissue
- · These Gauze Swab available in American fold and British fold
- · Gauze Swab available with X-Ray detectable thread

### **APPLICATIONS**

- It is an useful utility product in kit can be used for padding & protection
- · For swabbing and cleaning wounds
- · For dressing low exudates wounds



# Medical Roller Bandage

### **Absorbent Gauze Roll**

### SIZE:

 $90~\text{CM}~\text{X}~100~\text{MTR}~(1~\text{KG}~),~90~\text{CM}~\text{X}~91~\text{MTR}~(1.5~\text{KG}),\\ 90~\text{CM}~\text{X}~91~\text{MTR}~(2.3~\text{KG}),~90~\text{CM}~\text{X}~50~\text{MTR}$ 

### **CHARACTERISTICS**

- Made From Accomplished Cotton Yarn
- Absorption rate faster (less than 10 seconds)
- · Better tending capacity
- · Less allegiance to healing tissue
- These Absorbent Gauze Roll available in British Standard as well as United States Standard
- · Also available with X-ray detectable thread

- · High Utility & Life
- · Easy to Use
- These are uniquely designed to make them perfect for all types of application







# Medical Roller Bandage

### **M'PORE**

### Microporous Paper Tape

SIZE:

2.5 CM x 5 MTR, 5 CM x 5 MTR, 7.5 CM x 5 MTR, 10 CM x 5 MTR, 2.5 CM x 9.1 MTR, 5 CM x 9.1 MTR, 7.5 CM x 9.1 MTR, 10 CM x 9.1 MTR

### **CHARACTERISTICS**

- It is Pressure Sensitive Adhesive Tape
- · Hypoallergenic Adhesive, Breathable
- APPLICATIONS
- · Used In first aid to hold the dressing on to a wound



### Medica Rolled Bandage

### **CHARACTERISTICS**

- Made from superior cotton yarn AS PER SCH F (II) of I.P. / as per B.P.
- · Packed individually in plastic films

### **APPLICATIONS**

· For better dressing wounds



### **MEDITAPE**

### Adhesive Tape USP

### SI7F:

2.5~CM~x~5~MTR,~5~CM~x~5~MTR,~7.5~CM~x~5~MTR,~10~CM~x~5~MTR,~2.5~CM~x~8~Mtr~/~9~Mtr~/~10~Mtr,~5~CM~x~8~Mtr~/~9~Mtr~/~10~Mtr,~7.5~CM~x~8~Mtr~/~9~Mtr~/~10~Mtr,~7.5~CM~x~8~Mtr~/~9~Mtr~/~10~Mtr

### **CHARACTERISTICS**

- · Easy unrolling
- Strong adhesion
- · Easy & painless removal
- 20 % zinc oxide imparts it mild antiseptic, antiphlogistic, and antiinflammatory properties.
- Also provides the Hospital Packing like 8 Mtr, 9 Mtr, 10 Mtr.

### **APPLICATION**

· For dressing purpose







# Orthopaedic Products

### Gamjee Roll

### SIZE:

10 CM x 3 MTR, 15 CM x 3 MTR, 10 CM x 2 MTR, 15 CM x 2 MTR, 25 CM x 3 MTR, 25 CM x 2 MTR,

### **CHARACTERISTICS**

- Absorbent Cotton I.P. / B.P. Wrapped in Gauze Swab B.P.
- · Higher Absorption Rate

### **APPLICATIONS**

- For Swabbing, Cleaning Wound & Dressing Purpose
- · To Decrease Purulent Wound Drainage



### **MEDICREPE**

### Cotton Crepe Bandage B.P.

### SI7F

6 CM x 4 MTR (ST. LENGTH), 8 CM x 4 MTR (ST. LENGTH), 10 CM x 4 MTR (ST. LENGTH), 15 CM x 4 MTF H)

### **CHARACTERISTICS**

- · Better elasticity
- · Better aesthetic touch due to inter woven
- Better compliance
- Maintains uniform and controlled smooth pressure
- Reusable Washable
- APPLICATIONS
- · For strains, sprains, painful joints and after cast removal
- · Extreamly convenient as a pressure dressing for bums, skin grafts
- · General surgical, orthopadic, sport injuries and post operative conditions
- Very useful as a light compression bandage for muscular support



### **MEDIPOP**

### Plaster of Paris bandage B.P.

### SIZE:

5 CM x 2.7 MTR, 7.5 CM x 2.7 MTR, 10 CM x 2.7 MTR, 15 CM x 2.7 MTR

### **CHARACTERISTICS**

- · Abundant Creamy Texture
- · Outstanding Bandable
- · Heavy Duty Cast
- · Made as per British Pharmacopoiea
- APPLICATIONS
- · Immobilization Following Fractures
- · Orthopaedic Correction Treating Joints & Bone Conditions







# **Orthopaedic Products**

### **MEDIPLAST**

### Elastic Adhesive Bandage

### SIZE:

6 CM x 1 MTR (ST. LENGTH), 6 CM x 4/6 MTR (ST. LENGTH), 8 CM x 1 MTR (ST. LENGTH), 8 CM x 4/6 MTR (ST. LENGTH)

### **CHARACTERISTICS**

- · Better elasticity allows free muscle and joint movement
- Better aesthetic touch due to inter woven
- Porous adhesive permits skin breathing Strong adhesion
- Contains 20% Zinc Oxide which acts as a antiseptic, astringent and anti-inflammatory
- · Easy & painless removal



- · Provides the great support.
- · Immobilization and pressure where required.



### **MEDISOFT**

### Ortho Roll Cast Padding

### SI7F:

7.5 CM x 3 MTRS, 10 CM x 3 MTRS, 15 CM x 3 MTRS

### **CHARACTERISTICS**

- Made up of Non-Woven Viscose Material.
- Great Strength & Stretch Resistance Provides Good Support to Cast.
- Its work as a cushion between skin & cast

### **APPLICATION**

 Helps in Overcoming the Situation, when there is Oedema (Swelling) at the Fractured Site.



### **CANFIX**

### Cannula Fixator

### SI7F

MEDIUM (5 CM X 6 CM) (Individually Packed)

### **CHARACTERISTICS**

- · Simple to use, painless removal.
- · No irritation and no traces on the skin.
- · Allows air and humidity to pass freely.
- · Non-allergic adhesive with paper liner for easy removal.

### **APPLICATION**

• Fixes the cannula in the region of insertion till the procedure is completed.







# UAE \_

Atlantic World
Al Khaimah,UAE
Call: +971 52 540 1488
Address:Atlantic World FZC,Rakez,
Ras Al Khaimah-UAE.

# Atlantic World

INDIA =

New Delhi, India
Call: +91 6388 936 704
Address: D-130, Jaima Nagar.

NEW DELHI-INDIA.

# TURKEY G

Atlantic World Istanbul, Turkey Call: +905052713995 Address: Office-27,3443, Beyoglu/Isanbul, TURKEY.

# EGYPT ==

Atlantic World Cairo, Egypt Address : Office-545, Rehab City,New Cairo. EGYPT.

# NIGERIA .

Atlantic World & Global Investment Nigeria Limited-ABUJA-NIGERIA Call: +2349093221961

Address: MF-78,Top Hill,Abuja, F.C.T.

NIGERIA.



info@atlantic-world.com



www.atlanticworld.com